



# Kangmin Hygienic Material Development Co.,Ltd

Add.: Dingluan Industrial Park, Changyuan city, Henan Province, China

## Medical N95 respirator instructions

Type: Sterilized by EO. EO residues<10μg

Certificate: CE

Wear way: Ear loop

Material : Non-woven

Bacterial Filtration Efficiency >95%

Particle Filtration Efficiency: >95%

Usage: Medical protection,It is suitable for filtering particles in the air, blocking droplets, blood, body fluids and secretions under medical working environment

### Features

1. high quality material reach high fluid resistance.
2. Adjustable nose piece
- 3.Class I flammability.
4. high filtration efficiency.
5. EO residues<10μg.
- 6.very low breathing resistance.
- 7.Ultra-Sonic Welding way ensure good quality and strong mask wear

### How to use:

1. Open the package,
2. Keep the mask in the hand and mask with bridge is top,
3. Wear the ear-loop on the nose
4. Adjust the bridge until the mask fit well with nose and face.

### Packing Note

20pcs/inner box, 600pcs/master box by size: 45\*43.5\*41.5cm

Product: Medical protective mask

Type Designation: KMP-01




Figure 1 General view of mask  
(The information in this figure will be replaced by Copy of marking plate in final package)




Figure 2 General view of mask

Product: Medical protective mask

Type Designation: KMP-01




Figure 3 General view of mask



Figure 4 General view of mask

Product: Medical protective mask

Type Designation: KMP-01




Figure 5 General view of mask




Figure 6 General view of mask (4 ply)

END OF THE PHOTO DOCUMENTATION

